

MSL - CEILING MOUNT HEIGHT DETAILS FOR SERIES 1 & 3

To: Medical & Surgical Lighting


From:.....

Project:.....

Order #:.....

Provide measurements below on the ceiling type at site.

(B) Suspended Ceiling



NOTE:
 In installations where the combined distance between the "floor to false ceiling" and the "false ceiling to the fixing point" is greater than 3600mm, a separate mounting structure will be required to be installed by client.

See MSL "Cotton Reel Extension" option. Available as a separate PDF download and optional purchase.

Please contact MSL to discuss any queries you might have to achieve the best possible solution.

For Installations direct to a structural ceiling, please advise dimension **A**
 For Installations with a suspended ceiling, please advise dimensions at **B** and **C**

IMPORTANT:
 As all ceiling mounts are custom made to suit each installation, this form must be completed and returned to MSL prior to manufacture.

